

Mega Bed (Electric 3 Function)



A premium three-function electric hospital bed engineered for modern healthcare environments, combining robust cold-rolled steel construction with intuitive remote-controlled operation. Featuring three high-quality, noise-free Linak motors and a full range of adjustable positioning, it delivers effortless patient care and comfort — straight out of the box.

Product Features

- Three Linak motors, noise-free, remote control
- Auto-contour function for patient positioning
- Detachable HPL headboard/footboard, anti-fire resistant
- Two-section side rails, EN 60601-2-52 compliant
- Infusion holes and drainage hooks
- Height-adjustable IV pole included
- Emergency button included

CODE
TKMB3FBW



Product Specifications

- Backrest adjustment: 0–70° ($\pm 5^\circ$)
- Legrest adjustment: 0–40° ($\pm 5^\circ$)
- Height adjustment: 480–680mm (± 30 mm)
- Frame: Cold-rolled steel, electrostatic powder coating
- Castors: Four 125mm, central brake
- Mattress: 1980 x 890 x 80mm PU sponge



www.constar.com.au

