

Mega Bedside Cabinet



A practical and well-appointed bedside cabinet designed for hospital and healthcare settings. Constructed with an aluminium alloy pillar and HPL board, it combines durable materials with thoughtful storage solutions — keeping everything a patient needs within easy reach while maintaining a clean, professional appearance at bedside.

Product Features

- Smooth, silent sliding drawer
- Enclosed cabinet plus open storage compartment
- 3-sided tabletop railing prevents items falling
- Towel bars on both sides
- Noise-free castors, two with brake
- Available in 4 colours including Beech Wood

CODE
TKBSCBW



Product Specifications

- Material: Aluminium alloy pillar, HPL board, aluminium composite panel
- Storage: 1 drawer, 1 enclosed cabinet, 1 open compartment
- Castors: 4 x $\phi 60$ mm, 2 with brake
- Dimensions: L510 x W510 x H800mm



www.constar.com.au

