

Mega Overbed Table



A sleek and functional overbed table designed for patient comfort and clinical practicality. Built with an aluminium alloy lifting column and pneumatic height adjustment, it offers effortless one-handed operation alongside a 90° rotating HPL tabletop — making it easy to position for eating, reading, or working bedside.

Product Features

- Pneumatic height adjustment, one-handed operation
- HPL tabletop rotates 90° clockwise and counterclockwise
- Tabletop rail secures books, devices and accessories
- U-shaped base for easy cleaning and access
- 4 x noise-free castors, 2 with brake
- Available in Beech, Chocolate, Walnut

CODE
TKOBTBW



Product Specifications

- Material: Aluminium alloy column, cold-rolled steel base
- Height adjustment: 770–1,130mm
- Max load capacity: 10kg
- Castors: 4 x ϕ 50mm, 2 with brake
- Dimensions: L880 x W430 x H770–1,130mm



www.constar.com.au

