

Task ProCare Hot & Cold Therapy Sleeve S



Task ProCare Hot & Cold Therapy Sleeve is a reusable flexible gel sleeve designed to deliver targeted thermal therapy for ankle recovery and musculoskeletal support. Its dual-function design allows for both cold and heat application, with the flexible gel conforming naturally to the ankle for consistent, comfortable contact during treatment.

Product Information

- Reusable flexible gel sleeve
- Size: 23 x 12cm
- Fit: Ankle
- Dual-function: hot and cold therapy
- Conforms closely for consistent contact
- Microwave in 10-sec increments for heat
- Freeze 2+ hours; apply up to 15 minutes cold



Packaging

- 1 sleeve per box
- 28 boxes per carton
- Total pieces per carton - 28

Carton Size

- 38x27.5x27.5cm
- G.W 00kg
- N.W 00kg

Box Size

- 23x12.5x2.5cm
- N.W 00KG

Inner Box



Outer Carton



www.diagmed.com.au

